

Select "Scan" and use the iPad camera to select the barcode for Donor ID

Anesthesia Touch will recognize the barcode for Donor ID. Move the camera or product label to center the barcode.

Anesthesia Touch will confirm a correct scan and list the unscanned barcodes. Scan in the other barcodes in any order to complete the process.

You may also select “Close” (upper left) to cancel

After scanning the barcodes, set the time, rate, and comments for the blood product. Select “Save” when complete.

Blood products will appear on the record in the same way as any fluid. Select the product in the medication grid to modify or remove.

[illegible]