

Anesthesia Touch 4.5

Platform Requirements



iOS Device Requirements

Component	Requirement
Hardware version	All iPad versions supported, including Mini, Air, and Pro
Free Disk Space	10GB
Operating System	iOS 14 or Later
Wireless Networking	WPA 2 or WPA 2 Enterprise

Workstation

Hardware	Requirement
Processor	1 GHz (x86, x64 processor) or faster Recommended: Dual core processors 1.8GHz or faster
Memory	4 GB RAM or greater
Disk Space	Minimum: 100 GB or greater
Network Card	100Base-T; Wi-Fi also supported
Keyboard & Mouse	Standard Keyboard & Mouse. Optional medical grade and/or washable
Monitor	Touchscreen monitor is recommended for workstations in the OR, optimally 1280x1024 or larger.
Tablet Functionality / Form Factor	Microsoft Surface Pro 6 or Later
Printer	Either a network printer installed or an attached printer installed

Software	Requirement
Operating System	Microsoft Windows 10. 32-bit or 64-bit
Antivirus	Standard antivirus software compliant with facility's policy
Remote Desktop & Support software	SecureLink Quick Connect or Windows Remote Desktop via SecureLink
.Net Platform	Microsoft .net 4.0 client profile
Compatible Browsers	IE 10+/Edge, Chrome, Safari, Firefox (Latest or One Prior)
Adobe Reader	9 or above

Streaming Data Capture Requirements

Component	Requirement
Network connectivity	Ethernet Drop or Wi-Fi Network
Power	Power outlet in close proximity of anesthesia monitor or machine

Touch Connect Service

Touch Connect is a locally installed service that manages inbound and outbound interfaces for the Touch Platform as well as negotiating secure SSL connections to the Plexus Cloud. See the *Touch Connect Introduction and Requirements* for more information.

Touch Connect Requirements

Touch Connect has the following requirements and does not need a dedicated server. It can operate in parallel with any antivirus software installed.

Component	Requirement
Hard Drive Space	2 GB
Memory Recommended	4GB – 8GB
Processor	Dual Core
Network Connection	Minimum: 100 Base-T. Recommended: 1000 Base-T
Internet Connection	Broadband or LTE internet connection recommended
Touch Connect Ports	443 to the Plexus Cloud
SecureLink Gatekeeper	Ports 22 and 80 to the SecureLink server
Host OS Requirements	- Windows 7, 8, or 10 - Windows Server 2008 R2, 2012, 2016, 2019, or 2022

SecureLink Gatekeeper

For customer support, Plexus Technology Group uses SecureLink Gatekeeper for remote access. Plexus Technology Group will not require VPN accounts for access when SecureLink Gatekeeper is installed. Please see the SecureLink Gatekeeper documentation for more information.